# Chymotrypsin Aktivität KIT

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## 15. LITERATUR


# Chymotrypsin Activity KIT

*For the in vitro determination of Chymotrypsin activity in stool*

valid from January 26th, 2005

K 6990

## 1. INTENDED USE

The described photometric kit is intended for the activity determination of Chymotrypsin in stool. It is for *in vitro* diagnostic use only.

## 2. INTRODUCTION

### Indications

- Exocrine pancreas insufficiency•in act within chronic pancreatitis

### 3 MATERIAL SUPPLIED

Catalogue No	Content	Kit Components	Quantity
K 6990SUB	SUB	Substrate, lyophilized	20 x 2 ml
K 6990KA	CTRLA	Control A, lyophilized (range see specification)	4 x 500 µl
K 6990KB	CTRLB	Control B, lyophilized (range see specification)	4 x 500 µl

### 4 MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Precision pipettors and disposable tips to deliver 10-1000 µl
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- 10 mm one way cuvette
- Spectral photometer for measurement kinetics at 405 nm (Hg)
- Standard laboratory glass or plastic vials, cups, etc.
- Sample preparation system (e.g. Roche Catalog No. 745 804)

## 5. PREPARATION AND STORAGE OF REAGENTS

- The lyophilized **substrate** (SUB) is stable at **2-8 °C** until expiry date stated on the label. Substrate has to be reconstituted with **2000 µl** aqua dest. Allow to dissolve for 10 minutes and then mix thoroughly by gentle inversion to insure complete reconstitution. The reconstituted substrate is stable at 2-8 °C for 5 days and at 15- 25 °C for 2 days.
- The lyophilized **controls** (CTRL) are stable at **2-8 °C** until the expiry date stated on the label. Controls have to be reconstituted with **500 µl** aqua dest. Allow to dissolve for 10 minutes and then mix thoroughly by gentle inversion to insure complete reconstitution. Reconstituted controls are **not stable**.

## 6. SAMPLE PREPARATION

### *Extraction of the stool sample*

We recommend as **sample extraction buffer** the solvens buffer (Catalogue No. K 6990SOL). We recommend the following sample preparation:

We recommend the use of a stool sample preparation kit for dosing 100 mg of stool sample (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany; cat # 745804). The stool sample has to be suspended in 10 ml solvens buffer.

**Constant buffer volume: 10 ml**

**Constant dilution factor: 1:100**

In case the sample preparation kit is not used the stool sample weight has to be diluted with 100 times of solvens buffer (e.g. 120 mg stool sample + 120 ml of solvens buffer).

The stool suspension has to be mixed with the solvens (e.g. with a vortexer for 30 sec.).

Afterwards the suspension is transferred to a tube with a cap (e.g. Eppendorf tube) and centrifuged 5 minutes at 13000 g.

## *Dilution of samples*

### **Stool samples**

For analysis, pipette **100 µl** of the supernatant without dilution per cuvette.

**The stool suspension is not stable and can not be stored.** We recommend to weigh fresh sample amount for a new assay, if the analysis should be repeated.

The Chymotrypsin activity in stool is constant for 10 days at room temperature.

## **7. ASSAY PROCEDURE**

### *Principle of the test*

Succ-Ala-Ala-Pro-Phe-pNA + H<sub>2</sub>O  $\xrightarrow{\text{Chymotrypsin}}$  Succ-Ala-Ala-Pro-Phe + p-Nitroanilin

### *Test procedure*

Run the assay either at 25 °C or at 37 °C (see data sheet)

Add 100 µl of CTRL (control A, control B) or sample into the respective cuvette.

Add 2000 µl SUB (substrate) into each cuvette.

Determine the extinction after approximate 1 min at 405 nm and start the stop clock synchronically. Determine the extinction again after exactly 1, 2 and 3 min.

## 8. RESULTS

The mean is calculated from the extinction difference pro min ( $\Delta E/\text{min}$ )

Calculation:

The activity of the chymotrypsin from the sample is given in the table below or is calculated as in the following example described:

$$\text{U/g Stool} = 212 \times \Delta E_{405}/\text{min}$$

Values - table for measurement at Hg 405 nm

$\Delta E/\text{min}$	U/g	$\Delta E/\text{min}$	U/g	$\Delta E/\text{min}$	U/g
0,001	0,2	0,013	2,8	0,045	9,5
0,002	0,4	0,014	3,0	0,050	10,6
0,003	0,6	0,015	3,2	0,055	11,7
0,004	0,8	0,016	3,4	0,060	12,7
0,005	1,1	0,017	3,6	0,065	13,8
0,006	1,3	0,018	3,8	0,070	14,8
0,007	1,5	0,019	4,0	0,075	15,9
0,008	1,7	0,020	4,2	0,080	17,0
0,009	1,9	0,025	5,3	0,085	18,0
0,010	2,1	0,030	6,4	0,090	19,0
0,011	2,3	0,035	7,4	0,095	20,1
0,012	2,5	0,040	8,5	0,100	21,2

## 9. LIMITATIONS

The dilution limit is for a detection time of 3 min at a extinction difference of  $\Delta E/\text{min} = 0,100$ . For sample with higher activity mix 0.1 ml suspension with 0.4 ml of solvent buffer and re-assay the sample. The result has to be multiplied with 5.

## 10. QUALITY CONTROL

Immundiagnostik AG recommends the use of commercial control samples for internal quality control if available.

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Expected values*

#### Normal ranges

(1g stool corresponds to 1 ml 1ml)

	25°C	37°C
normal range	> 6 U/g	> 13,2 U/g
grey area	3 – 6 U/g	6,6 – 13,2 U/g
pathological range	< 3 U/g	< 6,6 U/g

Lit.: H.Goebbel et al. *Labormedizin*, Supplement 1 (1985) 8-10

We recommend each laboratory to establish its own norm concentration range.

## 11. PERFORMANCE CHARACTERISTICS

## 12. PRECAUTIONS

- For *in vitro* diagnostic use only.
- The quality control guidelines should be followed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

## 13. TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according to the enclosed manual.

## 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG together with a written complaint.

## 15. REFERENCES

21.03.2005 26012005Chymotrypsin.DOC