

A T C A T C G T A T C C A T C C C T A T C C A T C C C T
T C A T C C C T A T C G C T T C C A T C T G C T T C A
A T C T T C A T C C G T A T C T G C T C C A T C T G C T
C C A T T T C A T C C G T A T C
G C T T C C A T C C C T A T C G C
T C A T C C A T C C C T A T C T G C T T C C A T C T G C T



BIOTECON Diagnostics

foodproof[®]

foodproof[®]
D-Light

D-Light

Reagent D



For safer food - BIOTECON Diagnostics

Differentiation of DNA from dead and living cells

The **foodproof[®]** D-Light instrument is designed for the differentiation of DNA from dead and living cells. The reaction principle is a photoactivation using Reagent D to avoid false-positive PCR results. **Reagent D** contains a light sensitive substance which can only penetrate the cell membranes of dead cells. Exposure to visible light leads to covalent binding of this substance to DNA and prevents the DNA from being amplified via PCR. It works as a preliminary step of sample treatment before the DNA extraction. The stand-alone **foodproof[®]** D-Light is a compact benchtop solution that fits in every lab.

BIOTECON Diagnostics offers its customers complete workflow solutions, either via manual or automated sample preparation. The combination of the **foodproof[®]** D-Light and Reagent D, **foodproof[®]** DNA extraction kits and subsequent real-time PCR detection kits and all assay-required laboratory equipment represents a true “One-Stop-Shop” solution for the customer.

Fast: Up to 96 samples can be processed simultaneously

Convenient: Just insert plate and close lid, program starts automatically

Versatile: Insertion slide-in allows use of different microplates and adapter in SBS format e.g. deepwell plates

Safe: Secure sample workflow; no ice needed, as no development of heat; reduced risk of cross contamination

Economic: LED with very low energy consumption of just 9 W

Approved: Tested and certified for CE/FCC and product safety

Applications using **foodproof[®]** D-Light

- *Enterobacteriaceae*
- *Cronobacter*
- *Salmonella*
- *E. coli*
- *Shigella*
- *Vibrio*
- Yeast & Mold

Please contact BIOTECON Diagnostics for application details

A T C A T C G T A T C C A T C C C T A T C C A T C C C T
 T C A T C C C T A T C G C T T C C A T C T G C T T C A
 A T C T T C A T C C G T A T C T G C T C C A T C T G C T
 C C A T T T C A T C C G T A T C
 G C T T C C A T C C C T A T C G C
 T C A T C C A T C C C T A T C T G C T T C C A T C T G C T



BIOTECON Diagnostics

Workflow



**Enrichment/
Subcultivation**



**Reagent D
Treatment**



DNA Extraction



**Real-time PCR
Detection**

Specifications D - Light

| | |
|-------------------------------------|---|
| Dimensions (W x D x H) | 140 mm x 232 mm x 90 mm |
| Weight | 1.6 kg |
| Electrical data | 100-240 VAC, 50/60Hz, 1.0-0.5A, 15VDC, 0.8A |
| Power demand (standby) | 0.3 W |
| Power demand (while running) | 9.0 W |
| Working temperature range, humidity | 5 – 40 °C; 35 – 75 % rel. humidity |
| Storage temperature range, humidity | 0 – 65 °C; 35 – 75 % rel. humidity |
| Time set for dark and blue phase | 0 – 99 minutes for dark and |
| Acoustic signals | Signal after initialization and |
| Temperature control | Active control, shutdown at 50 °C |
| Voltage control | 14.2 to 15.8 V |
| Outage detection | LED outage |

Order Numbers

D 110 45 **foodproof®** D-Light
 A 500 02 Reagent D

Available through BIOTECON Diagnostics:

Robotic System

Applications and technical service

DNA Extraction & Detection Kits

For foodborne Pathogens, GMO, Allergen, Animal ID, etc.

Additional

All necessary reagents and consumables

BIOTECON Diagnostics GmbH

Hermannswerder 17
 14473 Potsdam
 Germany

Phone: +49 (0)331-2300-200
 Fax: +49 (0)331-2300-299

roboprep@bc-diagnostics.com
 www.bc-diagnostics.com