

STANDARD F COVID-19 Ag FIA

STANDARD F COVID-19 Ag FIA is the fluorescent immunoassay for the qualitative detection of specific antigens to SARS-CoV-2 present in human nasopharynx.

CE
Cat. No : 10COV30D

SD BIOSENSOR

Product Specification

- Test time** : Within 30 mins
- Specimen** : Nasopharyngeal swab or NP swabs collected into a viral transport medium
- Storage temperature** : 2-30°C/36-86°F
- Fluorescent Immunoassay (Europium)** : Higher sensitivity than rapid test (X 4times)
- Result Analyzer** : F2400, F200, F100
- Pack size** : 25 Tests/kit



Benefit

- Excellent sensitivity and specificity via FIA method
- Easy to use
- Fast results within 30 minutes
- More than 4 times higher sensitivity than RDTs
- Room temperature storage
- Ready-to-use reagents
- Automated platform with small POC analyzer

Performance Characteristics

Clinical evaluation

This clinical evaluation of STANDARD F COVID-19 Ag FIA was conducted with using 155 patients' NP swab specimens in viral transport media. The RT-PCR method was used as confirmatory test to reveal patients' status of SARS-CoV-2 infection.

1) Analysis by days after symptom onset

Days after symptom onset	Sensitivity	Positive	Negative	Total
≤ 7 days	89.36%	42	5	47
> 7 days	83.33%	5	1	6
asymptomatic	100.00%	2	0	2

2) Analysis by RT-PCR Ct value (RdRP gene basis)

Ct value	Sensitivity	Positive	Negative	Total
≤ 20	100.00%	11	0	11
20 < Ct ≤ 30	94.44%	34	2	36
> 30	50.00%	4	4	8

- 97.22% Positive agreement as Ct value of RT-PCR method was below 30


3) Clinical Sensitivity & Specificity

	RT-PCR (FDA EUA approved)		Total
	Positive	Negative	
Positive	49	4	53
Negative	6	96	102
Total	55	100	155
Clinical Sensitivity	89.09% (49/55), (95% CI : 77.75% - 95.89%)		
Clinical Specificity	96% (96/100), (95% CI : 90.07% - 98.90%)		

Test Procedure

'STANDARD TEST' Mode - For a single patient specimen


STANDARD F100, F200 and F2400 analyzer




1 Insert the test device to the test slot of the analyzer.



2 Apply 4 drops of extracted specimen.



3 Press the 'TEST START'



4 Result in 30 mins.

'READ ONLY' Mode - For multiple simultaneous tests

STANDARD F100 and F200 analyzer




1 Apply 4 drops of extracted specimen.



2 Incubate for 30 mins.



3 Set the analyzer to the 'Read Only' Mode.



4 Read the test result.

Specification of STANDARD F Analyzers



- F100** : Hand-held sized & Battery-powered POC analyzer
- F200** : Table top analyzer with maximized usability
- F2400** : High-throughput analyzer for mass screening test

STANDARD F Analyzer is a next-generation Fluorescent Immunoassay system. It measures multiple bio-markers within a single platform, and three different models are capable of covering various medical & laboratory settings.

Ordering Information

STANDARD F Analyzers

Cat. No.	Product	Unit	Weight	Carton size (W/D/H)
10FA10	F100	1 Unit	0.7 kg	105 x 135 x 100 mm
10FA20	F200	1 Unit	2.5 kg	200 x 240 x 205 mm
10FA24	F2400	1 Unit	20 kg	510 x 566 x 297 mm

STANDARD F COVID-19 Ag FIA

Cat. No.	Product	Storage temperature	Tests/kit	Kits/carton	Carton size (W/D/H)
10COV30D	STANDARD F COVID-19 Ag FIA	2-30°C/36-86°F	25	30	560 x 520 x 390 mm
10COVC10	STANDARD COVID-19 Ag Control	2-30°C/36-86°F	10	20	155 x 390 x 110 mm

STANDARD www.sdbiosensor.com

Head Office C-4&5Floor, 16, Deogyong-daero 155beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690, REPUBLIC OF KOREA
Tel +82-31-300-0400 | Fax +82-31-300-0499 | E-mail sales@sdbiosensor.com © 2020 SD BIOSENSOR. All rights reserved.

LEF-FCOVAG-EN_R.02 / 202007

SD BIOSENSOR