

## Immunolab Aflatoxin B1 Immunoaffinity Columns

Parameter	Cat. No.
Aflatoxin B1 Immunoaffinity Columns	AB1-A01

### Intended Use:

The Immunolab Aflatoxin B1 Immunoaffinity Columns are used for the purification and enrichment of Aflatoxin B1 from different matrices. The resulting eluate can be quantitatively analysed by various analytical techniques like the Immunolab Aflatoxin B1 ELISA or HPLC.

### Specifications:

Test principle:	Immunoaffinity chromatography
Format:	10 immunoaffinity columns per kit
Capacity:	> 200 µg
Recovery (spiked samples):	97%
Storage:	2-8°C