

Immunolab Aflatoxin (total) Immunoaffinity Columns

Parameter	Cat. No.
Aflatoxin (total) Immunoaffinity Columns	AFT-A01

Intended Use:

The Immunolab Aflatoxin (total) Immunoaffinity Columns are used for the purification and enrichment of Aflatoxin (total) from different matrices. The resulting eluate can be quantitatively analysed by various analytical techniques like the Immunolab Aflatoxin B1 and M1 ELISA or HPLC.

Specifications:

Test principle:	Immunoaffinity chromatography
Format:	10 immunoaffinity columns per kit
Capacity:	> 200 µg
Recovery (spiked samples):	97%
Storage:	2-8°C