

## Immunolab Zearalenone Immunoaffinity Columns

Parameter	Cat. No.
Zearalenone Immunoaffinity Columns	ZEA-A01

### Intended Use:

The Immunolab Zearalenone Immunoaffinity Columns are used for the purification and enrichment of Zearalenone from different matrices. The resulting eluate can be quantitatively analysed by various analytical techniques like the Immunolab Zearalenone ELISA or HPLC.

### Specifications:

<b>Test principle:</b>	Immunoaffinity chromatography
<b>Format:</b>	10 immunoaffinity columns per kit
<b>Capacity:</b>	> 2000 µg
<b>Recovery (spiked samples):</b>	97%
<b>Storage:</b>	2-8°C