

Immunolab Sesame ELISA kit

Parameter	Cat. No.
Sesame ELISA	SES-E01

Intended Use:

The Immunolab Sesame ELISA is a quick and sensitive method for the quantitative determination of sesame in foods.

Specifications:

Test principle:	Sandwich enzyme immunoassay	
Format:	Microtiter plate with 96 wells	
Incubation time:	60 min (20/20/20 min)	
Standards:	0; 2; 5; 15; 30 ppm	
Detection:	450 nm (HRP + TMB)	
Detection limit (LOD):	Soup	0.2 ppm
	Icecream	0.2 ppm
	Sausage	0.2 ppm
	Salad sauce	0.2 ppm
	Cracker	0.2 ppm
Cross-reactivities:	Oat	0.0003%
	Bean	0.0003%
Storage:	2–8°C	

Spike Solution

Parameter	Concentration	Cat. No.
Sesame	1000 ppm	SES-S01

Intended Use:

The Immunolab Sesame Spike Solution is intended for the execution of spiking experiments in combination with the Immunolab Sesame ELISA kit.