

Immunolab Beta-Lactoglobulin ELISA kit

Parameter	Cat. No.
Beta-Lactoglobulin ELISA	BLG-E01

Intended Use:

The Immunolab Beta-Lactoglobulin ELISA is a quick and sensitive method for the quantitative determination of Beta-Lactoglobulin in foods.

Specifications:

Test principle:	Sandwich enzyme immunoassay
Format:	Microtiter plate with 96 wells
Incubation time:	60 min (20/20/20 min)
Standards:	0; 10, 40; 100; 400 ppb
Detection:	450 nm (HRP + TMB)
Detection limit (LOD):	1.5 ppb
Cross-reactivities:	Ewe's (sheep's) milk < 0.2% Goat's milk < 0.002% Casein < 0.02%
Storage:	2–8°C

Spike Solution

Parameter	Concentration	Cat. No.
Beta-Lactoglobulin	10000 ppb	BLG-S01

Intended Use:

The Immunolab Beta-Lactoglobulin Spike Solution is intended for the execution of spiking experiments in combination with the Immunolab Beta-Lactoglobulin ELISA kit.