

## Immunolab Casein ELISA kit

Parameter	Cat. No.
Casein ELISA	CAS-E01

### Intended Use:

The Immunolab Casein ELISA is a quick and sensitive method for the quantitative determination of Casein in foods.

### Specifications:

<b>Test principle:</b>	Sandwich enzyme immunoassay
<b>Format:</b>	Microtiter plate with 96 wells
<b>Incubation time:</b>	60 min (20/20/20 min)
<b>Standards:</b>	0; 0.2; 0.6; 2; 6 ppm
<b>Detection:</b>	450 nm (HRP + TMB)
<b>Detection limit (LOD):</b>	0.04 ppm
<b>Cross-reactivities:</b>	Sheep's (ewe's) milk 1.2% Goat's milk 1.1%
<b>Storage:</b>	2–8°C

## Spike Solution

Parameter	Concentration	Cat. No.
Casein	1000 ppm	CAS-S01

### Intended Use:

The Immunolab Casein Spike Solution is intended for the execution of spiking experiments in combination with the Immunolab Casein ELISA kit.