

## Immunolab Deoxynivalenol ELISA kit

<

Parameter	Cat. No.
Deoxynivalenol ELISA	DON-E01

### Intended Use:

The Immunolab Deoxynivalenol ELISA is a sensitive method for the quantitative determination of Deoxynivalenol in grain, beer and soy.

### Specifications:

<b>Test principle:</b>	Competitive enzyme immunoassay
<b>Format:</b>	Microtiter plate with 96 tests
<b>Incubation time:</b>	60 min (30/15/15 min)
<b>Standards:</b>	0; 2; 8; 20; 40; 80 [ng/mL]
<b>Detection:</b>	450 nm (HRP + TMB)
<b>Detection limit:</b>	Grain 12.5 ng/kg Beer 5 ng/kg Soy 12.5 ng/kg
<b>Cross-reactivities:</b>	15-Acetyl-DON 0.4% 3-Acetyl-DON 800%
<b>Storage:</b>	2–8°C