

Immunolab Fish ELISA kit

Parameter	Cat. No.
Fish ELISA	FIS-E01

Intended Use:

The Immunolab Fish ELISA is a quick and sensitive method for the quantitative determination of fish in foods.

Specifications:

Test principle:	Sandwich enzyme immunoassay
Format:	Microtiter plate with 96 wells
Incubation time:	60 min (20/20/20 min)
Standards:	0; 4; 10; 40; 100 ppm
Detection:	450 nm (HRP + TMB)
Detection limit (LOD):	1,4 ppm
Storage:	2–8°C

Spike Solution

Parameter	Concentration	Cat. No.
Fish (Cod)	1000 ppm	FIS-S01

Intended Use:

The Immunolab Fish (Cod) Spike Solution is intended for the execution of spiking experiments in combination with the Immunolab Fish (Cod)ELISA kit.