

Immunolab Fumonisin ELISA kit

Parameter	Cat. NO.
Fumonisin ELISA	FUM-E01

Intended Use:

Der Immunolab Fumonisin ELISA is a sensitive method for the quantitative determination of Fumonisin in grain and maize.

Specifications:

Test principle:	Competitive enzyme immunoassay
Format:	Microtiter plate with 96 wells
Incubation time:	80 min (60/20 min)
Standards:	0; 0,5; 2; 5; 10; 25 [ng/mL]
Detection:	450 nm (HRP + TMB)
Detection limit:	Grain 22,5 µg/kg Maize 22.5 µg/kg
Storage:	2–8°C