

Immunolab Lysozyme ELISA kit

Parameter	Cat. No.
Lysozyme ELISA	LYS-E01

Intended Use:

The Immunolab Lysozyme ELISA is a quick and sensitive method for the quantitative determination of lysozyme in wine and other foods.

Specifications:

Test principle:	Sandwich enzyme immunoassay
Format:	Microtiter plate with 96 wells
Incubation time:	60 min (20/20/20 min)
Standards:	0; 25, 50; 100; 250 ppb
Detection:	450 nm (HRP + TMB)
Detection limit (LOD):	2 ppb
Storage:	2–8°C