

## Immunolab Ovalbumin ELISA kit

Parameter	Cat. No.
Ovalbumin ELISA	OVA-E01

### Intended Use:

The Immunolab Ovalbumin ELISA is a quick and sensitive method for the quantitative determination of ovalbumin in wine and other foods.

### Specifications:

<b>Test principle:</b>	Sandwich enzyme immunoassay
<b>Format:</b>	Microtiter plate with 96 wells
<b>Incubation time:</b>	60 min (20/20/20 min)
<b>Standards:</b>	0; 25, 100; 250; 500 ppb
<b>Detection:</b>	450 nm (HRP + TMB)
<b>Detection limit (LOD):</b>	4 ppb
<b>Storage:</b>	2–8°C

## Spike Solution

Parameter	Concentration	Cat. No.
Ovalbumin	10000 ppb	OVA-S01

### Intended Use:

The Immunolab Ovalbumin Spike Solution is intended for the execution of spiking experiments in combination with the Immunolab Ovalbumin ELISA kit.