

Immunolab Progesterone ELISA kit

Parameter	Cat. No.
Progesterone ELISA	PRO-E01

Intended Use:

The Immunolab Progesterone ELISA is a sensitive method for the quantitative determination of unphysiologically high Progesterone concentrations in blood plasma.

Specifications:

Test principle:	Competitive enzyme immunoassay
Format:	Microtiter plate with 96 wells
Incubation time:	140 min (60/60/20 min)
Standards:	0; 1; 5; 10; 50; 100 [ng/mL]
Detection:	450 nm (HRP + TMB)
Detection limit:	Blood plasma 0.1 µg/L
Cross-reactivities:	17a-Hydroxyprogesterone 5.0% Hydroxycortison 0.3% Pregnenolone 1.1% Pregnanolone 0.5%
Storage:	2–8°C