

Immunolab Soja (Soy) ELISA kit

Parameter	Cat. No.
Soja (Soy) ELISA	SOJ-E01

Intended Use:

The Immunolab Soja ELISA is a quick and sensitive method for the quantitative determination of soja (soy) in foods.

Specifications:

Test principle:	Sandwich enzyme immunoassay
Format:	Microtiter plate with 96 wells
Incubation time:	60 min (20/20/20 min)
Standards:	0; 40; 100; 400; 1000 ppb STI
Detection:	450 nm (HRP + TMB)
Detection limit (LOD):	16 ppb STI (soja trypsin inhibitor)
Cross-reactivities:	Sesam 0.0002%
Storage:	2–8°C

Spike Solution

Parameter	Concentration	Cat. No.
Soy Trypsin Inhibitor	10000 ppb	SOJ-S01

Intended Use:

The Immunolab Soy Trypsin Inhibitor Spike Solution is intended for the execution of spiking experiments in combination with the Immunolab Soja (Soy) ELISA kit.