

Immunolab T-2 Toxin RAPID ELISA kit

Parameter

Cat. NO.

T-2 Toxin RAPID ELISA

T2T-E03

Intended Use:

The Immunolab T-2 Toxin RAPID ELISA is a sensitive method for the quantitative determination of T-2 Toxin in cereals, beer, milk and meat.

Specifications:

Test principle: Competitive enzyme immunoassay

Format: Microtiter plate with 96 wells

Incubation time: 20 min (10/10 min)

Standards: 0; 17.5; 87.5; 350; 875; 1750 [ppb]

Detection: 450 nm (HRP + TMB)

Detection limit:	Wheat	10 ppb
	Rye	11 ppb
	Barley	16 ppb
	Corn	14 ppb
	Rice	11 ppb
	Meat	8 ppb
	Milk	5 ppb
	Beer	5 ppb