

## Immunolab Zearalenone ELISA kit

Parameter	Cat. No.
Zearalenone ELISA	ZEA-E01

### Intended Use:

The Immunolab Zearalenone ELISA is a sensitive method for the quantitative determination of Zearalenone in grain.

### Specifications:

<b>Test principle:</b>	Competitive enzyme immunoassay
<b>Format:</b>	Microtiter plate with 96 wells
<b>Incubation time:</b>	200 min (120/60/20 min)
<b>Standards:</b>	0; 0,5; 2; 20; 200 [ng/mL]
<b>Detection:</b>	450 nm (HRP + TMB)
<b>Detection limit:</b>	Grain 12.5 µg/kg
<b>Cross-reactivities:</b>	α-Zearalenol 11% β-Zearalenol 16% α-Zearalanol 3%
<b>Storage:</b>	2–8°C