

SSI® SALMONELLA "BIG FIVE" SEROTYPING KIT

STATENS
SERUM
INSTITUT




Supplied by:

**Oxford
BioSystems**



Tel: +44(0)1235 431390
sales@oxfordbiosystems.com
www.oxfordbiosystems.com



SSI Diagnostica
2 Herredsvejen
DK-3400 Hillerød
Denmark

T +45 4829 9178
F +45 4829 9179
@ ssidiagnostica@ssi.dk
w ssi.dk/ssidiagnostica

SALMONELLA "BIG FIVE" SEROTYPING KIT

Description

SSI Diagnostica has composed a kit, which contains all the *Salmonella* antisera needed to serotype five important public health related *Salmonella* serovars according to ISO/TR 6579-3 "Microbiology of food and animal feed - Horizontal method for the detection, enumeration and serotyping of *Salmonella*".

The 5 *Salmonella* serovars are:

- *Salmonella* Typhimurium
- *Salmonella* Enteritidis
- *Salmonella* Infantis
- *Salmonella* Virchow
- *Salmonella* Hadar

The kit includes all relevant O antisera, H antisera, and H phase inversion antisera. All the antisera are also sold individually.

Background

ISO/TR 6579-3 has been under development during 2013. It gives guidance on serotyping of the five most common *Salmonella* serovars which has been isolated from humans.

They can cause serious food poisoning. The increased occurrence of *Salmonella* infections within the last decade has accentuated the need for serotyping of *Salmonella* as a base for proper diagnosis, identification of sources of infection and food control.

Storage and Shelf Life

The antisera should be kept in a dark place at 2-8 °C. The stability of SSI Diagnostica antisera is minimum 2 years from date of shipment. The expiry dates are printed on the individual reagents.

Product Assortment

Product	Art.no.
SSI® Big Five <i>Salmonella</i> Serovars ISO/TR 6579-3 Kit.	89074

Or order the individual sera:

O antisera:	Art. no.
O:4	23839
O:8	40224
O:6,7,8	40221
O:9	40276
O:6 ₁	40273
O:7	40274
O:46	40288

H antisera:	Art. no.
H:q,s,t,p,u	22667
H:i	23850
H:r	40304
H:z ₁₀	40308
H:m	40332
H:x	40330
H:2	40340
H:5	40341

H phase inversion antisera:	Art. no.
SG 1	40363
SG 2	40364
SG 5	40367
SG 6	40368
H:r phase inversion	40381

Information and Ordering

SSI Diagnostica
2 Herredsvejen
3400 Hillerød
Denmark
T +45 4829 9178
F +45 4829 9179
@ ssidiagnostica@ssi.dk
w ssi.dk/ssidiagnosticas
🛒 shop.ssi.dk

