


Human, Mouse & Rat Singleplex ELISA Kits

Q-Plex Singleplex assays are achieved by placing two spots consisting of capture antibodies in a defined array to the bottom of each well of a 96-well plate. Along with the assay spots, there is one positive control spot for assuring proper assay procedure and for software overlay placement. The array design of the singleplex kit, allows the user to have intra-well replicates for additional statistical data. Using less than 50 µl of sample per well, up to 80 samples can be assayed for a single assay within 2.5 hours.



Singleplex Benefits

- ▶ Better sensitivity than traditional ELISAs.
- ▶ Intra-well precision from replicate spots provides greater confidence of results.
- ▶ Easy to use ELISA Protocol.

Human Singleplex ELISA Kits

Cat. No.	Product Name	Cat. No.	Product Name	Cat. No.	Product Name
RQS462149HU	Q-Plex™ Human Adiponectin Assay	RQS465649HU	Q-Plex™ Human IL-1ra Assay	RQS467149HU	Q-Plex™ Human MCP-2 Assay
RQS462249HU	Q-Plex™ Human Adipsin Assay	RQS465749HU	Q-Plex™ Human IL-2 Assay	RQS467249HU	Q-Plex™ Human MCP-3 Assay
RQS462349HU	Q-Plex™ Human Ang-2 Assay	RQS461049HU	Q-Plex™ Human IL-2 HS* Assay	RQS467349HU	Q-Plex™ Human MIG Assay
RQS462449HU	Q-Plex™ Human CD-14 Assay	RQS466249HU	Q-Plex™ Human IL-4 Assay	RQS467449HU	Q-Plex™ Human MIP-1α Assay
RQS462649HU	Q-Plex™ Human CD-26 Assay	RQS461249HU	Q-Plex™ Human IL-4 HS* Assay	RQS467549HU	Q-Plex™ Human MIP-1β Assay
RQS462549HU	Q-Plex™ Human CD-163 Assay	RQS466349HU	Q-Plex™ Human IL-5 Assay	RQS468149HU	Q-Plex™ Human MMP 1 Assay
RQS462749HU	Q-Plex™ Human Cortisol Assay	RQS461349HU	Q-Plex™ Human IL-5 HS* Assay	RQS467749HU	Q-Plex™ Human MMP 2 Assay
RQS462849HU	Q-Plex™ Human C-Peptide Assay	RQS466449HU	Q-Plex™ Human IL-6 Assay	RQS467849HU	Q-Plex™ Human MMP 3 Assay
RQS462949HU	Q-Plex™ Human CRP Assay	RQS461449HU	Q-Plex™ Human IL-6 HS* Assay	RQS467949HU	Q-Plex™ Human MMP 7 Assay
RQS463049HU	Q-Plex™ Human CTACK Assay	RQS466549HU	Q-Plex™ Human IL-6R Assay	RQS468049HU	Q-Plex™ Human MMP 9 Assay
RQS463149HU	Q-Plex™ Human CXCL5 Assay	RQS466649HU	Q-Plex™ Human IL-7 Assay	RQS467649HU	Q-Plex™ Human MMP 13 Assay
RQS463249HU	Q-Plex™ Human E1G Assay	RQS466749HU	Q-Plex™ Human IL-8 Assay	RQS468249HU	Q-Plex™ Human MPO Assay
RQS463349HU	Q-Plex™ Human E2 Assay	RQS464749HU	Q-Plex™ Human IL-10 Assay	RQS468349HU	Q-Plex™ Human P4 Assay
RQS463449HU	Q-Plex™ Human Eotaxin Assay	RQS460349HU	Q-Plex™ Human IL-10 HS* Assay	RQS468449HU	Q-Plex™ Human PAI-1 Assay
RQS469349HU	Q-Plex™ Human sFAS Assay	RQS464849HU	Q-Plex™ Human IL-12p40 Assay	RQS468549HU	Q-Plex™ Human PDG Assay
RQS463549HU	Q-Plex™ Human FGF Assay	RQS464949HU	Q-Plex™ Human IL-12p70 Assay	RQS468649HU	Q-Plex™ Human PDGF Assay
RQS463649HU	Q-Plex™ Human Fractalkine Assay	RQS460449HU	Q-Plex™ Human IL-12p70 HS* Assay	RQS468749HU	Q-Plex™ Human PF-4 Assay
RQS463749HU	Q-Plex™ Human FSH Assay	RQS465049HU	Q-Plex™ Human IL-13 Assay	RQS468849HU	Q-Plex™ Human P-Selectin Assay
RQS463849HU	Q-Plex™ Human G-CSF Assay	RQS460549HU	Q-Plex™ Human IL-13 HS* Assay	RQS468949HU	Q-Plex™ Human RANTES Assay
RQS463949HU	Q-Plex™ Human Ghrelin Assay	RQS465149HU	Q-Plex™ Human IL-15 Assay	RQS469049HU	Q-Plex™ Human RBP4 Assay
RQS464049HU	Q-Plex™ Human Glucagon Assay	RQS460649HU	Q-Plex™ Human IL-15 HS* Assay	RQS469149HU	Q-Plex™ Human Resistin Assay
RQS464149HU	Q-Plex™ Human GM-CSF Assay	RQS465249HU	Q-Plex™ Human IL-16 Assay	RQS469249HU	Q-Plex™ Human S100A9 Assay
RQS464249HU	Q-Plex™ Human GRO-α Assay	RQS465349HU	Q-Plex™ Human IL-17 Assay	RQS469449HU	Q-Plex™ Human TARC Assay
RQS464349HU	Q-Plex™ Human HCGb Assay	RQS460749HU	Q-Plex™ Human IL-17 HS* Assay	RQS469549HU	Q-Plex™ Human TIMP-1 Assay
RQS464449HU	Q-Plex™ Human HGF Assay	RQS465849HU	Q-Plex™ Human IL-21 Assay	RQS469649HU	Q-Plex™ Human TIMP-2 Assay
RQS464549HU	Q-Plex™ Human I-309 Assay	RQS465949HU	Q-Plex™ Human IL-22 Assay	RQS469749HU	Q-Plex™ Human TNF-α Assay
RQS464649HU	Q-Plex™ Human IFN-γ Assay	RQS466049HU	Q-Plex™ Human IL-23 Assay	RQS461549HU	Q-Plex™ Human TNF-α HS* Assay
RQS460249HU	Q-Plex™ Human IFN-γ HS* Assay	RQS461149HU	Q-Plex™ Human IL-23 HS* Assay	RQS469849HU	Q-Plex™ Human TNF-β Assay
RQS465449HU	Q-Plex™ Human IL-1α Assay	RQS466149HU	Q-Plex™ Human IL-27 Assay	RQS461649HU	Q-Plex™ Human TNF-β HS* Assay
RQS460849HU	Q-Plex™ Human IL-1α HS* Assay	RQS466849HU	Q-Plex™ Human IP-10 Assay	RQS469949HU	Q-Plex™ Human TNF-RI Assay
RQS465549HU	Q-Plex™ Human IL-1β Assay	RQS466949HU	Q-Plex™ Human Leptin Assay	RQS460049HU	Q-Plex™ Human TNF-RII Assay
RQS460949HU	Q-Plex™ Human IL-1β HS* Assay	RQS467049HU	Q-Plex™ Human MCP-1 Assay	RQS460149HU	Q-Plex™ Human VEGF Assay

* High sensitivity detection

Mouse Singleplex ELISA Kits

Cat. No.	Product Name	Cat. No.	Product Name	Cat. No.	Product Name
RQS462149MS	Q-Plex™ Mouse Eotaxin Assay	RQS462949MS	Q-Plex™ Mouse IL-5 Assay	RQS463749MS	Q-Plex™ Mouse MDC Assay
RQS462249MS	Q-Plex™ Mouse GM-CSF Assay	RQS463049MS	Q-Plex™ Mouse IL-6 Assay	RQS463849MS	Q-Plex™ Mouse MIP-1 α Assay
RQS462349MS	Q-Plex™ Mouse IFN- γ Assay	RQS463149MS	Q-Plex™ Mouse IL-10 Assay	RQS463949MS	Q-Plex™ Mouse MIP-2 Assay
RQS462449MS	Q-Plex™ Mouse IL-1 α Assay	RQS463249MS	Q-Plex™ Mouse IL-12 Assay	RQS464049MS	Q-Plex™ Mouse RANTES Assay
RQS462549MS	Q-Plex™ Mouse IL-1 β Assay	RQS463349MS	Q-Plex™ Mouse IL-13 Assay	RQS464149MS	Q-Plex™ Mouse TARC Assay
RQS462649MS	Q-Plex™ Mouse IL-2 Assay	RQS463449MS	Q-Plex™ Mouse IL-17 Assay	RQS464249MS	Q-Plex™ Mouse TCA-3 Assay
RQS462749MS	Q-Plex™ Mouse IL-3 Assay	RQS463549MS	Q-Plex™ Mouse KC Assay	RQS464349MS	Q-Plex™ Mouse TNF- α Assay
RQS462849MS	Q-Plex™ Mouse IL-4 Assay	RQS463649MS	Q-Plex™ Mouse MCP-1 Assay		

Rat Singleplex ELISA Kits

Cat. No.	Product Name	Cat. No.	Product Name
RQS462149RT	Q-Plex™ Rat IFN- γ Assay	RQS462649RT	Q-Plex™ Rat IL-6 Assay
RQS462249RT	Q-Plex™ Rat IL-1 α Assay	RQS462749RT	Q-Plex™ Rat IL-10 Assay
RQS462349RT	Q-Plex™ Rat IL-1 β Assay	RQS462849RT	Q-Plex™ Rat IL-12p70 Assay
RQS462449RT	Q-Plex™ Rat IL-2 Assay	RQS462949RT	Q-Plex™ Rat TNF- α Assay
RQS462549RT	Q-Plex™ Rat IL-4 Assay		

› Manufactured by



Quansys Biosciences
365 North 500 West, Logan, UT 84321

› www.quansysbio.com

› Distributed by



BioVendor - Laboratorni medicina a.s.

Karasek 1767/1
621 00 Brno, Czech Republic
Phone: +420 549 124 185
Fax: +420 549 211 460
E-mail: info@biovendor.com

› www.biovendor.com

BioVendor GmbH

Otto-Hahn-Straße 16
34123 Kassel, Germany
Phone: +49 6221 4339 100
Fax: +49 6221 4339 111
E-mail: infoEU@biovendor.com

BioVendor GesmbH

Nußdorfer Straße 20/10
1090 Vienna, Austria
Phone: +43 1 890 9025
Fax: +43 1 890 9025-15
E-mail: infoAustria@biovendor.com



Oxford Biosystems Ltd

115J Olympic Avenue, Milton Park, Oxfordshire
OX14 4SA, United Kingdom
Phone: 01235 431390, E-mail: sales@oxfordbiosystems.com

› www.oxfordbiosystems.com